

**Measurements of complex dielectric constants of samples
in millimeter wave bands using Martin-Puplett type
Fourier Transform Spectrometer with high sensitive
Millimeter-wave bolometers**

Instructor: Makoto Hattori (hattori@astr.tohoku.ac.jp)

GPPU Experimental Point (GEP): 4

Goal of Study

Learn following things. 1. Fundamentals and operation procedures of the Martin-Puplett type Fourier Transform Spectrometer (MP-FTS) [1]. 2. Methods for measuring sensitivity and fundamentals of semiconductor bolometer detector system [2]. 3. How to make 0.28K environment with He3 sorption refrigerator. 4. Treatments of Turbomolecular Pump system and Mechanical Booster Pump system. 5. Treatments of wet dewar. 6. Data acquisition (DAQ) systems: read out circuit with operational amplifier and Labview. 7. How to measure the Frequency dependence of the complex dielectric constants of samples in millimeter wavebands using the FTS. 8. Writing data analysis codes.

Contents


MP-FTS is a broad band absolute spectrometer which has played important role in wide field of astronomy, cosmology and calibrating detector systems in the wave bands longer than far infrared wave bands. The FIRAS instrument mounted on COBE satellite [3] is one of the most famous applications of MP-FTS. It measured spectrum of the cosmic microwave background (CMB) radiation and showed that it follows the perfect black body spectrum with temperature of 2.725K. This result has confirmed that our universe has been evolved following the description of the standard big band theory. A photo of the MP-FTS which is used in this experiment is shown in right panel. MP-FTS does not lose value and is still active in observational astronomy and cosmology fields.



In this course, the application of the MP-FTS to the measurements of the frequency dependence of complex dielectric constants of samples against millimeter wave at room temperature and liquid Nitrogen temperature is experienced. Millimeter wave high sensitive bolometer with Nitrogen Trans Doped Germanium (NTD-Ge) semiconductor thermistor is used as detector. A photo of inside of detector cryostat, which is wet dewar, is shown in left panel. Scenery of the measurement is shown in right panel. An example of the sample measurement auto-correlation fringe is shown



in left panel. An example of the sample measurement auto-correlation fringe is shown




in left panel. An example of the

Obtained frequency dependence of the transmittance of a sample against millimeter waves is shown in right panel. Beat signals found in the fringe and sinusoidal pattern found in the transmittance are caused by Fabry-Perot interference effect.



By fitting the transmittance, refractive index and absorption coefficient, that is real and imaginary part of dielectric constant, are obtained.



Textbook and References

- [1] **MP-FTS:** Martin,D.H. and Pulett,E., Infrared Phys., 10, 105 (1969)
- [2] **Detector:** G.H.Rieke, Detection of light, Cambridge Univ. Press, 1994
- [3] **Application example for cosmology:** J.C.Mathor et al., ApJ, 420, 439 (1994)
- [4] **Optics:** M.Born and E.Wolf, Principles of Optics, Sections 1 and 7 (1974)
- [5] **Example:** S.Adachi, M.Hattori, et al., Review of Scientific Instruments, V.91, id.016103 (2020)
- [6] **Cross-Dragon:** C.Dragone, The Bell system technical journal, vol.57., no.7. , 2663-2683(1978)

Progress Schedule

- ✧ **Day 1**
Study fundamentals and operation procedures of Semiconductor bolometers including read our circuit, He3 sorption refrigerator and MP-FTS with observing devices.
Close the cryostat.
- ✧ **Day 2**
Start pumping the cryostat with turbo pump. Precool the cryostat using liquid Nitrogen (LN2) after it reaches high enough vacuum degree.
While waiting that the cryostat reaches high enough vacuum degree, studies principles of extracting frequency dependence of complex dielectric constants from the data obtained by FTS and sensitivities of bolometer.
- ✧ **Day 3**
Stop the precooling processes and get out LN2 from the He layer of the cryostat.
Pour liquid He4 in the cryostat.
By decreasing pressure inside of He4 layer with the Mechanical Booster pump, lower the temperature down to cooler than the boiling temperature of He3.
By manipulating the temperature controller, cool down the detector mount to 0.28K.
Acquisite data to estimate detector sensitivities.
Perform measurements to get transmittance of sample against millimeter wave at room temperature and LN2 temperature, that is 77K.
Close the measurement system.
- ✧ **Days 4**
Analyze data obtained previous day.
Submit a report of the experiment and give a brief presentation on the experiment.
Discuss future possibilities of FTS. Fundamental idea and present status of the state of art technology 'on chip FTS' is introduced.

Other Details

| | |
|---------------------------|--|
| Course Period | 2021 Summer |
| Place | Science C building N402 |
| Number of Students | 1—3 |
| Evaluation method | The evaluation method will be based on experiments (60 %), presentation (30 %), and discussion (10 %). Single report is expected to be submitted from all participant but contribution of each person must be clarified. |

In Addition

M2 students who are nominated as GPPU students, are able to take this lecture.

In case F2F training is not possible due to COVID, change the theme as follows.

1. Optical simulation using simulator CODE-V. Attendance is requested to design optical system based cross-Dragonian telescope. The cross-Dragon system is the optical system which is able to compensate the instrumental polarization imprinted by the reflection of light by the primary mirror and is widely applied in CMB polarization experiments.

2. Develop the Markov Chain Monte Carlo (MCMC) based data analysis code to deduce the refractive index and the absorption coefficient of the material in millimeter wavebands. The refractive index and absorption coefficient of some material are deduced by applying the code to the data obtained by FTS system in our laboratory.

Schedule: Attendance must receive the CODE-V license key dongle and windows note book machine prior to the training course. First day is assigned to the lecture on cross-Dragon telescope and manipulation of CODE-V. Second day is assigned to the lecture of our FTS measurement system and MCMC. Third day is assigned to present the report of MCMC code for simple example. All the results are reported on final day. The first three days lectures are continuously performed. After two weeks interval, the final day's program is performed. During this interval, attendance must work individually in home.